



**PRODUCT SPECIFICATIONS**



## RELIABLE Physician's Chair Scales for Patients with Special Health Care Needs



Footrest lifts out of the way for safe seating and convenient transport of scale while the lift-away arms allow patients side or front access. Oversized wheels facilitate movement of scale. Complete with heavy-duty understructure, these scales come fully assembled and ready to use. Black vinyl upholstery is complemented by the metal support structure.

### MECHANICAL CHAIR SCALE

Mechanical chair scales are available in pounds or pounds/kilo models. Die-cast, easy-to-read beams are readable from front and back.

**Capacity - 475:** 400 lb x 4 oz  
**4751:** 400 lb x 4 oz/180 kg x 100 g

### DIGITAL CHAIR SCALE

Digital chair scales are available in pound or kilo models. Easy-to-read digital display rapidly shows weight to an accuracy of .2 lb (.1 kg). Battery power enhances portability.

**Capacity - 6475:** 400 lb x .2 lb  
**6475KG:** 180 kg x .1 kg  
**Display:** .7" high-contrast LCD  
**Power:** 6 "C" size Alkaline batteries  
(not included)  
**Option:** 115 VAC power adapter



Cardinal Scale reserves the right to improve, modify, or enhance features and specifications without prior notice.



**ONE DIGITAL SCALES** .com

Call Us Toll-Free

**800-620-7811**

[www.onedigitalscales.com](http://www.onedigitalscales.com)