




**BE ADAM**

## MDW Medical Scales



Speed, Performance, Value - Perfect Balance

- **Easy to move** With built in wheels and sturdy mechanical operations, with no need for any external power supply or batteries, this mechanical scale can be used anywhere from the locker room to the doctor's offices.
- **Safe non-slip platform** A durable molded non-slip platform will help protect patients when standing on the unit. Platform is removable for easy cleaning.
- **Reaching new heights** Measure the height of the person while on the unit to reduce the time and space needed to record vital information.
- **Plus these standard features:** Simple set-up and operation, wheels for mobility, clear and easy-to-read graduations.



The MDW is a completely portable weighing solution for medical professionals. A durable mechanical medical scale with an adjustable height rod which is ideal for doctors offices, locker rooms or home use. With its strong metal construction and wheels for mobility, this mechanical medical scale is built for moving. Ideal for doctor's offices, locker rooms, and home use the MDW's adjustable height rod and non-slip platform make taking measurements safe and quick.



### Technical Specifications:

<b>Model</b>	<b>MDW 160M</b>
<b>Capacity</b>	160Kg / 350Lb
<b>Readability</b>	0.1kg / 0.25Lb (4oz)
<b>Repeatability (s.d)</b>	0.2kg / 0.5Lb
<b>Linearity</b>	0.2kg / 0.5Lb
<b>Height range</b>	75-200cm / 29.5"-78.7"
<b>Division of measurement</b>	0.5cm / 0.2"
<b>Platform Size</b>	375 x 275mm / 14.8" x 10.8"

### General Specifications:

<b>Weighing units</b>	Kg / Lb
<b>Applications</b>	Weighing
<b>Features</b>	Removable anti-slip platform Adjustable zero point Wheels to allow easy relocation of unit Height rod
<b>Tare range</b>	Full capacity by subtraction
<b>Power supply</b>	No external power required (mechanical operation)
<b>Display</b>	Graduated markings
<b>Housing</b>	Metal construction
<b>Net weight</b>	15Kg / 33.1Lb
<b>Overall dimensions (w x d x h)</b>	530 x 275 x 1485mm / 20.9" x 10.8" x 58.5"

Adam Equipment has been providing professionals with quality balances and scales for over 30 years. The company is certified to ISO 9001-2000.

Adam Equipment manufactures and distributes a complete selection of mass measurement products for the laboratory, medical, educational, industrial and retail markets. All Adam balances and scales are backed by Adam assurance warranties with responsive and reliable service before, during and after the sale.

**ONE DIGITAL SCALES**.com

Your First Choice  
for Industrial, Home and Specialty Scales

**800-620-7811**

www.onedigitalscales.com

