


COMPACT SCALE

NEW

Full digital calibration &
A standard battery pack
for complete portability.
Calibration weight provided

A choice of units of measure:
Grams, Ounces,
Pennyweights, Troy Ounces



MP SERIES


MP SERIES

Model	Capacity	Readability	Pan Size ϕ
MP 200*	200 gm	0.01 gm	120mm / 4.7inch
MP 300*	300 gm	0.01 gm	
MP 600	600 gm	0.1 gm	155 x 155mm / 6.1 x 6.1inch
MP 1000	1000 gm	0.1 gm	
MP 1500	1500 gm	0.1 gm	
MP 2000H	2000 gm	0.1 gm	
MP 2000	2000 gm	1 gm	
MP 5000	5000 gm	1 gm	

* with Acrylic windshield



Model	MP 200*	MP 300*	MP 600	MP 1000	MP 1500	MP 2000H	MP 2000	MP 5000
Capacity (gm)	200	300	600	1000	1500	2000	2000	5000
Accuracy (gm)	0.01		0.1			1		
Weighing Units	Gram, oz, dwt, ozt							
Tareing Range	100%							
Calibration	Auto							
Stability Time	1.5 Seconds							
Auto Power Off	Yes							
Temperature	0 - 40°C (32 - 104°F)							
Pan Size (mm / in)	120 / 4.7		155 x 155 / 6.1 x 6.1					
Display	4 digit LCD Display							
Weight	475 gms		600 gm					
Package	Standard carton : 24.5 x 16.5 x 5.2 (cm ³) 12 Units in one box : 50 x 17.5 x 32 (cm ³)							
Power Source	6 X AA Size batteries or AC/DC (Not Included) Adapter 9V100mA							




APPLICATION :

- * Widely used in laboratory, medical & clinical use
- * Weighing of gold, silver & precious stones
- * Weighing of herbal medicines & grains
- * Letter / Postal mail weighing
- * Diet Mixing

FEATURES :

- * Most Economical.
- * Compact & Space saving.
- * Auto calibration & Auto zero tracking.
- * Backlight LCD Display for easy viewing.
- * Low Battery indication.
- * Weighing units gm, ct, dwt, oz, ozt available.



WORLDWIDE MARKETING:

ONE DIGITAL SCALES.com

Call Us Toll-Free

800-620-7811

www.onedigitalscales.com

Email :