



MS20 Series Medical / Veterinary Scales

Features

The MS20 scale is particularly suited to the medical or veterinary industry. With a capacity of 44 lb (20 kg) and readability of 1/2 oz (10 g), the MS20 will accurately weigh infants or small animals.

The scale has a "HOLD" feature that will capture an unstable weight and keep it on the display for recording after the weight is taken off.

The MS20 is offered in two models and can operate on battery (single 9VDC) or AC (adapter included). The MS20 comes with a plastic weighing tray and the MS20S comes with a stainless steel weighing tray.

Specifications

Capacity-

<u>lb</u>	<u>kg</u>
44 lb x 1/2 oz	20 kg x 10 g

Weight Units- lb / kg

Power Supply- 9V alkaline battery
or 9V 100mA AC adapter

Tray- Stainless steel or plastic

Inside Tray Dimensions-

MS20 (plastic)	20 1/2" x 10" x 2 3/8"
MS20S (stainless steel)	20 1/2" x 11 1/2" x 2 1/4"

Display- 1.0" LCD display with 4 active digits

Operating Temperature- 41° to 95°F (5° to 35°C)

Function Keys- ON/TARE/OFF, UNIT, HOLD

Approvals- Full FCC compliance
UL/CUL power adapter



MS20



MS20S



Your First Choice
for Industrial, Home and Specialty Scales

800-620-7811

www.onedigitalscales.com