

# PRECISION BALANCE

## Model VB-301 / VB-302

### Features

- Auto shut off – selectable (VB-302 only)
- Auto zero tracking
- Auto calibration
- Auto backlight (optional)
- Low battery indication
- Large LCD display (6 digits 15mm high)
- Stainless steel weighing pan
- kg/g/lb/oz (VB-301)
- g/oz/lb/ tl.T/fl.H/ct/gn/pcs (VB-302)
- 1/5000 ~ 1/50000 division available

### Specifications

Model No	Capacity	Graduation	Platter size	Windshield
VB-302-300	300g	0.01g	∅ 100mm	Yes
VB-302-500	500g	0.01g	∅ 100mm	Yes
VB-302-1500	1500g	0.1g	140x165mm	No
VB-302-3000	3000g	0.1g	140x165mm	No
VB-301-3000	3000g	1g	140x165mm	No
VB-301-6000	6000g	1g	140x165mm	No
VB-301-10000	10000g	2g	140x165mm	No
Shipping Weight	1kg			
Operating Temp.	0 ~ 40 °C			
Power Source	6 x AA Batteries (not included) or AC/DC Adaptor 9-12V / 100mA			

\* Specifications are subject to change without prior notice



VB-302



VB-301